

3-oxo-5-alpha-steroid 4-dehydrogenase 1

Catalog No: #AP87711



Package Size: #AP87711-1 50ug #AP87711-2 100ug #AP87711-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | 3-oxo-5-alpha-steroid 4-dehydrogenase 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | SR type 1, Steroid 5-alpha-reductase 1, S5AR 1 |
| Accession No. | Uniprot ID:Q68FF9 |
| Uniprot | Q68FF9 |
| GeneID | 78925; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Srd5a1

Note: This product is for in vitro research use only