

Neuronal-specific septin-3

Catalog No: #AP87713



Package Size: #AP87713-1 50ug #AP87713-2 100ug #AP87713-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuronal-specific septin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	3-Sep
Accession No.	Uniprot ID:Q9Z1S5
Uniprot	Q9Z1S5
GeneID	24050;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: 43711

Note: This product is for in vitro research use only