

Solute carrier family 12 member 2

Catalog No: #AP87714



Package Size: #AP87714-1 50ug #AP87714-2 100ug #AP87714-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 12 member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Basolateral Na-K-Cl symporter, Bumetanide-sensitive sodium-(potassium)-chloride cotransporter 1, NKCC1
Accession No.	Uniprot ID:P55011
Uniprot	P55011
GeneID	6558;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC12A2

Note: This product is for in vitro research use only