

Alpha-synuclein

Catalog No: #AP87715



Package Size: #AP87715-1 50ug #AP87715-2 100ug #AP87715-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-synuclein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Accession No.	Uniprot ID:Q3I5G7
Uniprot	Q3I5G7
GeneID	641350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SNCA

Note: This product is for in vitro research use only