

SAM pointed domain-containing Ets transcription factor

Catalog No: #AP87716

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87716-1 50ug #AP87716-2 100ug #AP87716-3 1mg

Description

| | |
|--------------------|--|
| Product Name | SAM pointed domain-containing Ets transcription factor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Prostate epithelium-specific Ets transcription factor,Prostate-specific Ets,Prostate-derived Ets factor,PDEF,PSE |
| Accession No. | Uniprot ID:O95238 |
| Uniprot | O95238 |
| GeneID | 25803; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SPDEF

Note: This product is for in vitro research use only