

Suppressor of cytokine signaling 1

Catalog No: #AP87719



Package Size: #AP87719-1 50ug #AP87719-2 100ug #AP87719-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Suppressor of cytokine signaling 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SOCS-1
Accession No.	Uniprot ID:Q9QX78
Uniprot	Q9QX78
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Socs1

Note: This product is for in vitro research use only