

Suppressor of cytokine signaling 2

Catalog No: #AP87721



Package Size: #AP87721-1 50ug #AP87721-2 100ug #AP87721-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Suppressor of cytokine signaling 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SOCS-2,Cytokine-inducible SH2 protein 2,CIS-2,STAT-induced STAT inhibitor 2,SSI-2,CIS2,SSI2,STATI2
Accession No.	Uniprot ID:O14508
Uniprot	O14508
GeneID	8835;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SOCS2

Note: This product is for in vitro research use only