

Calcipressin-2

Catalog No: #AP87724



Package Size: #AP87724-1 50ug #AP87724-2 100ug #AP87724-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcipressin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Down syndrome candidate region 1-like 1, Myocyte-enriched calcineurin-interacting protein 2, MCIP2, Regulator of calcineurin 2, Thyroid hormone-responsive protein ZAKI-4, DSCR1L1, ZAKI4
Accession No.	Uniprot ID:Q14206
Uniprot	Q14206
GeneID	10231;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RCAN2

Note: This product is for in vitro research use only