

Regulatory-associated protein of mTOR

Catalog No: #AP87727



Package Size: #AP87727-1 50ug #AP87727-2 100ug #AP87727-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Regulatory-associated protein of mTOR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Raptor,p150 target of rapamycin (TOR)-scaffold protein,Raptor
Accession No.	Uniprot ID:Q8K4Q0
Uniprot	Q8K4Q0
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rptor

Note: This product is for in vitro research use only