

SH3 domain-binding glutamic acid-rich-like protein 3

Catalog No: #AP87728



Package Size: #AP87728-1 50ug #AP87728-2 100ug #AP87728-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SH3 domain-binding glutamic acid-rich-like protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q91VW3
Uniprot	Q91VW3
GeneID	73723;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sh3bgrl3

Note: This product is for in vitro research use only