

# Protein S100-A1

Catalog No: #AP87730



Package Size: #AP87730-1 50ug #AP87730-2 100ug #AP87730-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein S100-A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	S-100 protein alpha chain,S-100 protein subunit alpha,S100 calcium-binding protein A1
Accession No.	Uniprot ID:P35467
Uniprot	P35467
GeneID	295214;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:S100a1

Note: This product is for in vitro research use only