

## Protein S100-A12

Catalog No: #AP87732



Package Size: #AP87732-1 50ug #AP87732-2 100ug #AP87732-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein S100-A12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Calgranulin-C,CAGC,Extracellular newly identified RAGE-binding protein,EN-RAGE,S100 calcium-binding protein A12
Accession No.	Uniprot ID:P80310
Uniprot	P80310
GeneID	100301483;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:S100A12

Note: This product is for in vitro research use only