

Selenocysteine insertion sequence-binding protein 2-like

Catalog No: #AP87736

Package Size: #AP87736-1 50ug #AP87736-2 100ug #AP87736-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Selenocysteine insertion sequence-binding protein 2-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SECIS-binding protein 2-like,KIAA0256
Accession No.	Uniprot ID:Q93073
Uniprot	Q93073
GeneID	9728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SECISBP2L

Note: This product is for in vitro research use only