

Secreted frizzled-related protein 2

Catalog No: #AP87742



Package Size: #AP87742-1 50ug #AP87742-2 100ug #AP87742-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secreted frizzled-related protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	sFRP-2, Protein SDF5, Secreted apoptosis-related protein 1, SARP-1, Sarp1, Sdf5
Accession No.	Uniprot ID:P97299
Uniprot	P97299
GeneID	20319;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sfrp2

Note: This product is for in vitro research use only