

Adenosylhomocysteinase

Catalog No: #AP87750



Package Size: #AP87750-1 50ug #AP87750-2 100ug #AP87750-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosylhomocysteinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AdoHcyase,S-adenosyl-L-homocysteine hydrolase,SAHH
Accession No.	Uniprot ID:P23526
Uniprot	P23526
GeneID	191;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AHCY

Note: This product is for in vitro research use only