

Regulator of G-protein signaling 19

Catalog No: #AP87751



Package Size: #AP87751-1 50ug #AP87751-2 100ug #AP87751-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Regulator of G-protein signaling 19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	RGS19
Accession No.	Uniprot ID:Q08DC7
Uniprot	Q08DC7
GeneID	539036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RGS19

Note: This product is for in vitro research use only