

Secretogranin-3

Catalog No: #AP87752



Package Size: #AP87752-1 50ug #AP87752-2 100ug #AP87752-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secretogranin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	1B1075, Secretogranin III, SgIII
Accession No.	Uniprot ID:P47868
Uniprot	P47868
GeneID	116635;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Scg3

Note: This product is for in vitro research use only