

## Stromal cell-derived factor 1

Catalog No: #AP87754



Package Size: #AP87754-1 50ug #AP87754-2 100ug #AP87754-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Stromal cell-derived factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SDF-1,12-O-tetradecanoylphorbol 13-acetate repressed protein 1,TPAR1,C-X-C motif chemokine 12,Pre-B cell growth-stimulating factor,PBSF,Thymic lymphoma cell-stimulating factor,TLSF,Sdf1
Accession No.	Uniprot ID:P40224
Uniprot	P40224
GeneID	20315;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cxcl12

Note: This product is for in vitro research use only