

## Alkaline phosphatase, tissue-nonspecific isozyme

Catalog No: #AP87756



Package Size: #AP87756-1 50ug #AP87756-2 100ug #AP87756-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alkaline phosphatase, tissue-nonspecific isozyme
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AP-TNAP,Alkaline phosphatase 2,Alkaline phosphatase liver/bone/kidney isozyme,Akp-2,Akp2
Accession No.	Uniprot ID:P09242
Uniprot	P09242
GeneID	11647;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Alpl

Note: This product is for in vitro research use only