

Sodium-dependent serotonin transporter

Catalog No: #AP87759



Package Size: #AP87759-1 50ug #AP87759-2 100ug #AP87759-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Sodium-dependent serotonin transporter |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | SERT,5HT transporter,5HTT,Solute carrier family 6 member 4 |
| Accession No. | Uniprot ID:P31652 |
| Uniprot | P31652 |
| GeneID | 25553; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc6a4

Note: This product is for in vitro research use only