

ZFHX4 Antibody

Catalog No: #44448

Package Size: #44448-1 50ul #44448-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZFHX4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ZFHX4 Antibody detects endogenous levels of total ZFHX4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZFHX4
Target Name	ZFHX4
Other Names	Zinc finger homeobox protein 4 Zinc finger homeodomain protein 4 ZFH-4
Accession No.	Swiss-Prot#: Q86UP3NCBI Gene ID:79776
Uniprot	Q86UP3
GeneID	79776;
Calculated MW	394 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only