

Retinoid isomerohydrolase

Catalog No: #AP87760



Package Size: #AP87760-1 50ug #AP87760-2 100ug #AP87760-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoid isomerohydrolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	All-trans-retinyl-palmitate hydrolase,Retinal pigment epithelium-specific 65 kDa protein,Retinol isomerase
Accession No.	Uniprot ID:Q9TVB8
Uniprot	Q9TVB8
GeneID	403803;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RPE65

Note: This product is for in vitro research use only