

Retinoic acid receptor RXR-gamma

Catalog No: #AP87761



Package Size: #AP87761-1 50ug #AP87761-2 100ug #AP87761-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor RXR-gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Nuclear receptor subfamily 2 group B member 3, Retinoid X receptor gamma, Nr2b3
Accession No.	Uniprot ID:Q5BJR8
Uniprot	Q5BJR8
GeneID	83574;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rxrg

Note: This product is for in vitro research use only