

## Synaptosomal-associated protein 25

Catalog No: #AP87768



Package Size: #AP87768-1 50ug #AP87768-2 100ug #AP87768-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Synaptosomal-associated protein 25
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SNAP-25, Synaptosomal-associated 25 kDa protein
Accession No.	Uniprot ID:Q17QQ3
Uniprot	Q17QQ3
GeneID	540853;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SNAP25

Note: This product is for in vitro research use only