

Selenoprotein S

Catalog No: #AP87771



Package Size: #AP87771-1 50ug #AP87771-2 100ug #AP87771-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Selenoprotein S
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SelS,Sg2,Vimp
Accession No.	Uniprot ID:Q8VHV8
Uniprot	Q8VHV8
GeneID	286900;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Selenos

Note: This product is for in vitro research use only