

Serpin A12

Catalog No: #AP87775



Package Size: #AP87775-1 50ug #AP87775-2 100ug #AP87775-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serpin A12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Visceral adipose tissue-derived serine protease inhibitor,Vaspin,Visceral adipose-specific serpin
Accession No.	Uniprot ID:Q7TMF5
Uniprot	Q7TMF5
GeneID	68054;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpina12

Note: This product is for in vitro research use only