

Serine/threonine-protein kinase SIK2

Catalog No: #AP87779



Package Size: #AP87779-1 50ug #AP87779-2 100ug #AP87779-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase SIK2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Salt-inducible kinase 2, SIK-2, Serine/threonine-protein kinase SNF1-like kinase 2, Snf1lk2
Accession No.	Uniprot ID:Q8CFH6
Uniprot	Q8CFH6
GeneID	235344;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SIK2

Note: This product is for in vitro research use only