

SUMO-activating enzyme subunit 1

Catalog No: #AP87780



Package Size: #AP87780-1 50ug #AP87780-2 100ug #AP87780-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SUMO-activating enzyme subunit 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ubiquitin-like 1-activating enzyme E1A,Aos1,Sua1,Ubl1a1,Uble1a
Accession No.	Uniprot ID:Q9R1T2
Uniprot	Q9R1T2
GeneID	56459;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sae1

Note: This product is for in vitro research use only