

## SUMO-activating enzyme subunit 1

Catalog No: #AP87781



Package Size: #AP87781-1 50ug #AP87781-2 100ug #AP87781-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SUMO-activating enzyme subunit 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ubiquitin-like 1-activating enzyme E1A,Uble1a
Accession No.	Uniprot ID:Q6AXQ0
Uniprot	Q6AXQ0
GeneID	308384;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sae1

Note: This product is for in vitro research use only