

# Sodium-dependent dopamine transporter

Catalog No: #AP87782



Package Size: #AP87782-1 50ug #AP87782-2 100ug #AP87782-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sodium-dependent dopamine transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DA transporter,Solute carrier family 6 member 3
Accession No.	Uniprot ID:P23977
Uniprot	P23977
GeneID	24898;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Slc6a3

Note: This product is for in vitro research use only