

Prosaposin

Catalog No: #AP87787



Package Size: #AP87787-1 50ug #AP87787-2 100ug #AP87787-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prosaposin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Sulfated glycoprotein 1,SGP-1,Sgp1
Accession No.	Uniprot ID:P10960
Uniprot	P10960
GeneID	25524;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Psap

Note: This product is for in vitro research use only