

Slit homolog 2 protein

Catalog No: #AP87788



Package Size: #AP87788-1 50ug #AP87788-2 100ug #AP87788-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Slit homolog 2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Slit-2
Accession No.	Uniprot ID:Q9WVC1
Uniprot	Q9WVC1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slit2

Note: This product is for in vitro research use only