

Serine--pyruvate aminotransferase, mitochondrial

Catalog No: #AP87791



Package Size: #AP87791-1 50ug #AP87791-2 100ug #AP87791-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine--pyruvate aminotransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SPT,Alanine--glyoxylate aminotransferase,AGT,Agxt1
Accession No.	Uniprot ID:O35423
Uniprot	O35423
GeneID	11611;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Agxt

Note: This product is for in vitro research use only