

Spectrin alpha chain, brain

Catalog No: #AP87793



Package Size: #AP87793-1 50ug #AP87793-2 100ug #AP87793-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Spectrin alpha chain, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-II spectrin,Fodrin alpha chain,Spectrin, non-erythroid alpha chain,Spna2,Spta2
Accession No.	Uniprot ID:P16086
Uniprot	P16086
GeneID	64159;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sptan1

Note: This product is for in vitro research use only