

# Regulator of G-protein signaling 19

Catalog No: #AP87794



Package Size: #AP87794-1 50ug #AP87794-2 100ug #AP87794-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Regulator of G-protein signaling 19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	RGS19,G-alpha-interacting protein,GAIP,GAIP,GNAI3IP
Accession No.	Uniprot ID:P49795
Uniprot	P49795
GeneID	10287;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:RGS19

Note: This product is for in vitro research use only