

Regulator of G-protein signaling 9

Catalog No: #AP87795



Package Size: #AP87795-1 50ug #AP87795-2 100ug #AP87795-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Regulator of G-protein signaling 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	RGS9
Accession No.	Uniprot ID:O46469
Uniprot	O46469
GeneID	281453;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RGS9

Note: This product is for in vitro research use only