

## Mothers against decapentaplegic homolog 2

Catalog No: #AP87797



Package Size: #AP87797-1 50ug #AP87797-2 100ug #AP87797-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mothers against decapentaplegic homolog 2   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Ms  |
| Other Names        | MAD homolog 2, Mad-related protein 2, mMad2, SMAD family member 2, SMAD 2, Madh2, Madr2 |
| Accession No.      | Uniprot ID:Q62432   |
| Uniprot            | Q62432  |
| GeneID             | 17126;  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                  |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Smad2

Note: This product is for in vitro research use only