

Runt-related transcription factor 2

Catalog No: #AP87804



Package Size: #AP87804-1 50ug #AP87804-2 100ug #AP87804-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Runt-related transcription factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Acute myeloid leukemia 3 protein,Core-binding factor subunit alpha-1,CBF-alpha-1,Oncogene AML-3,Osteoblast-specific transcription factor 2,OSF-2,Polyomavirus enhancer-binding protein 2 alpha A subunit,PEA2-alpha A,SL3-3 enhancer factor 1 alpha A subunit,SL3/AKV core-binding factor alpha A subunit,Aml3,Cbfa1,Osf2,Pebp2a
Accession No.	Uniprot ID:Q08775
Uniprot	Q08775
GeneID	12393;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Runx2

Note: This product is for in vitro research use only