

NAD-dependent protein deacylase sirtuin-5, mitochondrial

Catalog No: #AP87807

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87807-1 50ug #AP87807-2 100ug #AP87807-3 1mg

Description

Product Name	NAD-dependent protein deacylase sirtuin-5, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q3ZBQ0
Uniprot	Q3ZBQ0
GeneID	507347;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SIRT5

Note: This product is for in vitro research use only