

Sarcosine dehydrogenase, mitochondrial

Catalog No: #AP87809



Package Size: #AP87809-1 50ug #AP87809-2 100ug #AP87809-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sarcosine dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SarDH,BPR-2,DMGDHL1
Accession No.	Uniprot ID:Q9UL12
Uniprot	Q9UL12
GeneID	1757;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SARDH

Note: This product is for in vitro research use only