

Selenoprotein S

Catalog No: #AP87812



Package Size: #AP87812-1 50ug #AP87812-2 100ug #AP87812-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Selenoprotein S
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SelS, Minor histocompatibility antigen H47, H47, Vimp
Accession No.	Uniprot ID:Q9BCZ4
Uniprot	Q9BCZ4
GeneID	109815;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Selenos

Note: This product is for in vitro research use only