

S-phase kinase-associated protein 2

Catalog No: #AP87813



Package Size: #AP87813-1 50ug #AP87813-2 100ug #AP87813-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-phase kinase-associated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cyclin-A/CDK2-associated protein p45,F-box protein Skp2,F-box/LRR-repeat protein 1,p45skp2,FBXL1
Accession No.	Uniprot ID:Q13309
Uniprot	Q13309
GeneID	6502;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SKP2

Note: This product is for in vitro research use only