

Heterogeneous nuclear ribonucleoprotein A1

Catalog No: #AP87815



Package Size: #AP87815-1 50ug #AP87815-2 100ug #AP87815-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	hnRNP A1, Helix-destabilizing protein, Single-strand RNA-binding protein, Unwinding protein 1, UP1, hnRNP core protein A1, HNRPA1
Accession No.	Uniprot ID:P09867
Uniprot	P09867
GeneID	505093;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HNRNPA1

Note: This product is for in vitro research use only