

Sentan

Catalog No: #AP87816



Package Size: #AP87816-1 50ug #AP87816-2 100ug #AP87816-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sentan
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein S100-A1-like,S100 calcium-binding protein A1-like,S100A1L
Accession No.	Uniprot ID:A6NMZ2
Uniprot	A6NMZ2
GeneID	132203;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SNTN

Note: This product is for in vitro research use only