

Sodium-dependent dopamine transporter

Catalog No: #AP87822



Package Size: #AP87822-1 50ug #AP87822-2 100ug #AP87822-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium-dependent dopamine transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DA transporter,Solute carrier family 6 member 3,DAT1
Accession No.	Uniprot ID:Q01959
Uniprot	Q01959
GeneID	6531;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC6A3

Note: This product is for in vitro research use only