

Secretogranin-2

Catalog No: #AP87823



Package Size: #AP87823-1 50ug #AP87823-2 100ug #AP87823-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secretogranin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chromogranin-C, Secretogranin II, SgII, CHGC
Accession No.	Uniprot ID:P13521
Uniprot	P13521
GeneID	7857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SCG2

Note: This product is for in vitro research use only