

Heterogeneous nuclear ribonucleoproteins A2/B1

Catalog No: #AP87825



Package Size: #AP87825-1 50ug #AP87825-2 100ug #AP87825-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoproteins A2/B1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	hnRNP A2/B1,Hnrpa2b1
Accession No.	Uniprot ID:O88569
Uniprot	O88569
GeneID	53379;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hnrpa2b1

Note: This product is for in vitro research use only