

## DNA-directed RNA polymerase III subunit RPC9

Catalog No: #AP87826



Package Size: #AP87826-1 50ug #AP87826-2 100ug #AP87826-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA-directed RNA polymerase III subunit RPC9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	RNA polymerase III subunit C9, Calcitonin gene-related peptide-receptor component protein, CGRP-RCP
Accession No.	Uniprot ID:Q8VHM6
Uniprot	Q8VHM6
GeneID	114205;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Crp

Note: This product is for in vitro research use only