

## Sperm surface protein Sp17

Catalog No: #AP87834



Package Size: #AP87834-1 50ug #AP87834-2 100ug #AP87834-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sperm surface protein Sp17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cancer/testis antigen 22,CT22,Sp17-1,Sperm autoantigenic protein 17,Sperm protein 17,SP17
Accession No.	Uniprot ID:Q15506
Uniprot	Q15506
GeneID	53340;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SPA17

Note: This product is for in vitro research use only