

Palmitoyl-protein thioesterase 1

Catalog No: #AP87837



Package Size: #AP87837-1 50ug #AP87837-2 100ug #AP87837-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Palmitoyl-protein thioesterase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PPT-1,Palmitoyl-protein hydrolase 1,Cln1,Ppt
Accession No.	Uniprot ID:P45479
Uniprot	P45479
GeneID	29411;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppt1

Note: This product is for in vitro research use only